

RAMPE ALTA PRESSIONE

High pressure manifolds



	CODICE CODE	VIE WAYS	TIPO TYPE	P. MAX	ROTAZIONE ROTATION	CONNESSIONI CONNECTIONS	
RAMPE - 2 RUBINETTI Manifolds - 2 Stopcocks	M063009	2	RIGHT, OFF	600 PSI 41 BAR	180°	3 x L.L.F. ROTATING L.L.M.	100
	M063009/1	2	RIGHT, ON	600 PSI 41 BAR	180°	3 x L.L.F. ROTATING L.L.M.	100
	M063193	2	LEFT, OFF	600 PSI 41 BAR	180°	3 x L.L.F. ROTATING L.L.M.	100
	M063193/1	2	LEFT, ON	600 PSI 41 BAR	180°	3 x L.L.F. ROTATING L.L.M.	100
RAMPE - 3 RUBINETTI Manifolds - 3 Stopcocks	M062988	3	RIGHT, OFF	600 PSI 41 BAR	180°	4 x L.L.F. ROTATING L.L.M.	100
	M062988/1	3	RIGHT, ON	600 PSI 41 BAR	180°	4 x L.L.F. ROTATING L.L.M.	100
	M063866	3	LEFT, OFF	600 PSI 41 BAR	180°	4 x L.L.F. ROTATING L.L.M.	100
	M063866/1	3	LEFT, ON	600 PSI 41 BAR	180°	4 x L.L.F. ROTATING L.L.M.	100

Legenda: - L.L.M. = Luer Lock Maschio
- L.L.F. = Luer Lock Femmina

Legend: - **L.L.M. = Male Luer Lock**
- **L.L.F. = Female Luer Lock**